



Raised Toilet Seat Fitting Instructions

Reviewed: June 2025 Next review: June 2026 Version 1

About this leaflet

This leaflet provides instructions on how to safely fit your raised toilet seat at home.

What is the toilet seat for?

A raised toilet seat increases the height of your toilet, making it easier and safer to sit down and stand up, especially if you have reduced mobility or are recovering from surgery.

Fitting Instructions



Lift the toilet lid and seat.



Place the raised toilet seat directly onto the toilet bowl, as shown in the picture.



Push the front of the seat firmly against the front of the toilet bowl, ensuring it fits as closely as possible.



Tighten both side screws to secure the seat snugly against the bowl.

Check that the seat is securely fitted before each use, and regularly check for stability.

Contact Information

If you have any difficulty fitting the seat, please contact:

Occupational Therapy Department – 0141 951 5161

All of our publications are available in different languages, larger print, braille (English only), audio tape or another format of your choice.

我們所有的印刷品均有不同語言版本、大字體版本、盲文(僅有英文)、錄音 帶版本或你想要的另外形式供選擇。

كافة مطبو عاتنا مناحة بلغات مختلفة و بالأحرف الطباعية الكبيرة و بطريقة بريل الخاصة بالمكفوفين (باللغة الإنكليزية فقط) و على شريط كاسيت سمعي أو بصيغة بديلة حسب خيارك.

Tha gach sgrìobhainn againn rim faotainn ann an diofar chànanan, clò nas motha, Braille (Beurla a-mhàin), teip claistinn no riochd eile a tha sibh airson a thaghadh.

हमारे सब प्रकाशन अनेक भाषाओं, बड़े अक्षरों की छपाई, बेल (केवल अंग्रेज़ी), सुनने वाली कसेट या आपकी पसंदनुसार किसी अन्य फॉरमेट (आस्प) में भी उपलब्ध हैं। 我们所有的印刷品均有不同语言版本、大字体版本、盲文(仅有英文)、录音 带版本或你想要的另外形式供选择。

ਸਾਡੇ ਸਾਰੇ ਪਰਚੇ ਅਤੇ ਕਿਤਾਬਚੇ ਵਗ਼ੈਰਾ ਵੱਖ ਵੱਖ ਭਾਸ਼ਾਵਾਂ ਵਿਚ, ਵੱਡੇ ਅੱਖਰਾਂ ਅਤੇ ਬ੍ਰੇਲ (ਸਿਰਫ਼ ਅੰਗਰੇਜ਼ੀ) ਵਿਚ, ਅੱਡੀਓ ਟੇਪ 'ਤੇ ਜਾਂ ਤੁਹਾਡੀ ਮਰਜ਼ੀ ਅਨੁਸਾਰ ਕਿਸੇ ਹੋਰ ਰੂਪ ਵਿਚ ਵੀ ਮਿਲ ਸਕਦੇ ਹਨ।

ہاری تمام مطبوعات مختلف زبانوں، بوے حروف کی چھپائی، بریل (صرف اگریزی)، سنے والی کسٹ یا آپ کی پیند کے مطابق کسی دیگر صورت (فارمیٹ) میں بھی وستیاب ہیں۔



